

产品名称: FBP3 Rabbit Polyclonal Antibody
产品货号: APRab10861



产品概述 (Summary)

产品名称 (Production Name)	FBP3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FUBP3
别名 (Alternative Names)	FUBP3; FBP3; Far upstream element-binding protein 3; FUSE-binding protein 3
基因 ID (Gene ID)	8939.0
蛋白 ID (SwissProt ID)	Q96I24.The antiserum was produced against synthesized peptide derived from human FUBP3. AA range:201-250

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	60kDa

研究背景 (Background)


产品名称: FBP3 Rabbit Polyclonal Antibody
产品货号: APRab10861




function: May interact with single-stranded DNA from the far-upstream element (FUSE). May activate gene expression., PTM: Phosphorylated upon DNA damage, probably by ATM or ATR., similarity: Contains 4 KH domains., tissue specificity: Detected in a number of cell lines., function: May interact with single-stranded DNA from the far-upstream element (FUSE). May activate gene expression., PTM: Phosphorylated upon DNA damage, probably by ATM or ATR., similarity: Contains 4 KH domains., tissue specificity: Detected in a number of cell lines.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from COLO and HepG2 cells, using FBP3 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using FBP3 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA) .

产品名称: FBP3 Rabbit Polyclonal Antibody
产品货号: APRab10861



Immunohistochemical analysis of paraffin-embedded Human brain. Antibody was diluted at 1:100 (4°,overnight) . High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negative contrl (right) obtaned from antibody was pre-absorbed by immunogen peptide.

注意事项 (Note)

For research use only .