

产品名称: Fgl2 Rabbit Polyclonal Antibody
产品货号: APRab10953



产品概述 (Summary)

产品名称 (Production Name)	Fgl2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FGL2
别名 (Alternative Names)	FGL2; Fibroleukin; Fibrinogen-like protein 2; pT49
基因 ID (Gene ID)	10875.0
蛋白 ID (SwissProt ID)	Q14314.The antiserum was produced against synthesized peptide derived from human Fgl2. AA range:38-87

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	55kDa

研究背景 (Background)


产品名称: Fgl2 Rabbit Polyclonal Antibody
产品货号: APRab10953




fibrinogen like 2(FGL2) Homo sapiens The protein encoded by this gene is a secreted protein that is similar to the beta- and gamma-chains of fibrinogen. The carboxyl-terminus of the encoded protein consists of the fibrinogen-related domains (FRED). The encoded protein forms a tetrameric complex which is stabilized by interchain disulfide bonds. This protein may play a role in physiologic functions at mucosal sites. [provided by RefSeq, Jul 2008],function:May play a role in physiologic lymphocyte functions at mucosal sites.,similarity:Contains 1 fibrinogen C-terminal domain.,subunit:Homotetramer; disulfide-linked.,tissue specificity:Constitutively expressed in cytotoxic T-cells.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysate from 293 cells treated with insulin, using Fgl2 antibody.



Western Blot analysis of various cells using Fgl2 Polyclonal Antibody

注意事项 (Note)

For research use only .