

产品名称: FoxD4/D4L Rabbit Polyclonal Antibody
产品货号: AP Rab11082



产品概述 (Summary)

产品名称 (Production Name)	FoxD4/D4L Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FOXD4 FOXD4; FKHL9; FOXD4A; FREAC5; Forkhead box protein D4; Forkhead-related protein FKHL9; Forkhead-related transcription factor 5; FREAC-5; Myeloid factor-alpha; FOXD4L4; FOXD4B; Forkhead box protein D4-like 4; FOXD4-like 4; Forkhead box prote
别名 (Alternative Names)	
基因 ID (Gene ID)	653404/2298/100036519/349334/653427 Q12950/Q8WXT5/Q5VV16/Q6VB85/Q6VB84/Q3SYB3. The antiserum was
蛋白 ID (SwissProt ID)	produced against synthesized peptide derived from human FOXD4/L2/L3/L4/L5/L6. AA range:281-330

产品应用 (Application)

产品名称: FoxD4/D4L Rabbit Polyclonal Antibody
产品货号: AP Rab11082




稀释比 (Dilution Ratio) WB 1:500-1:2000, ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight) 46kDa

研究背景 (Background)

This gene encodes a member of the forkhead/winged helix-box (FOX) family of transcription factors. FOX transcription factors play critical roles in the regulation of multiple processes including metabolism, cell proliferation and gene expression during ontogenesis. Mutations in this gene are associated with a complex phenotype consisting of dilated cardiomyopathy, obsessive-compulsive disorders, and suicidality. [provided by RefSeq, Mar 2012], similarity: Contains 1 fork-head DNA-binding domain.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from Jurkat cells, using FOXD4/L2/L3/L4/L5/L6 Antibody. The lane on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .