

产品名称: FoxN2 Rabbit Polyclonal Antibody
产品货号: APRab11096



产品概述 (Summary)

产品名称 (Production Name)	FoxN2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FOXN2
别名 (Alternative Names)	FOXN2; HTLF; Forkhead box protein N2; Human T-cell leukemia virus enhancer factor
基因 ID (Gene ID)	3344.0
蛋白 ID (SwissProt ID)	P32314.The antiserum was produced against synthesized peptide derived from human FOXN2. AA range:341-390

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	36kDa

产品名称: FoxN2 Rabbit Polyclonal Antibody
产品货号: AP Rab11096




研究背景 (Background)


This gene encodes a forkhead domain binding protein and may function in the transcriptional regulation of the human T-cell leukemia virus long terminal repeat. [provided by RefSeq, Jul 2008],function: Binds to the purine-rich region in HTLV-I LTR.,similarity: Contains 1 fork-head DNA-binding domain.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from HeLa, 293, and HUVEC cells, using FOXN2 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from K562 cells using FOXN2 antibody.

产品名称: FoxN2 Rabbit Polyclonal Antibody
产品货号: APRab11096



Western Blot analysis of HT-29/NIH-3T3 using FoxN2 Polyclonal Antibody.

注意事项 (Note)

For research use only .