

产品名称: FPRL2 Rabbit Polyclonal Antibody
产品货号: APRab11125



产品概述 (Summary)

产品名称 (Production Name)	FPRL2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FPR3
别名 (Alternative Names)	FPR3; FPRH1; FPRL2; N-formyl peptide receptor 3; FMLP-related receptor II; FMLP-R-II; Formyl peptide receptor-like 2
基因 ID (Gene ID)	2359.0
蛋白 ID (SwissProt ID)	P25089.The antiserum was produced against synthesized peptide derived from human FPRL2. AA range:304-353

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC/IF 1:200-1:1000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	

产品名称: FPRL2 Rabbit Polyclonal Antibody
产品货号: APRab11125



研究背景 (Background)

function:Low affinity receptor for N-formyl-methionyl peptides, which are powerful neutrophils chemotactic factors.

Binding of FMLP to the receptor causes activation of neutrophils. This response is mediated via a G-protein that activates a phosphatidylinositol-calcium second messenger system.,similarity:Belongs to the G-protein coupled receptor 1


family.,function:Low affinity receptor for N-formyl-methionyl peptides, which are powerful neutrophils chemotactic factors.

Binding of FMLP to the receptor causes activation of neutrophils. This response is mediated via a G-protein that activates a phosphatidylinositol-calcium second messenger system.,similarity:Belongs to the G-protein coupled receptor 1 family.,

研究领域 (Research Area)

Neuroactive ligand-receptor interaction;

图片 (Image Data)



Immunofluorescence analysis of LOVO cells, using FPRL2 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .