

产品名称: G3BP2 Rabbit Polyclonal Antibody
产品货号: APRab11213



产品概述 (Summary)

产品名称 (Production Name)	G3BP2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	G3BP2
别名 (Alternative Names)	G3BP2; KIAA0660; Ras GTPase-activating protein-binding protein 2; G3BP-2; GAP SH3 domain-binding protein 2
基因 ID (Gene ID)	9908.0
蛋白 ID (SwissProt ID)	Q9UN86.The antiserum was produced against synthesized peptide derived from human G3BP2. AA range:433-482

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	55kDa

产品名称: G3BP2 Rabbit Polyclonal Antibody
产品货号: APRab11213




研究背景 (Background)


function: Probable scaffold protein that may be involved in mRNA transport .,PTM: Arg-457 and Arg-468 are dimethylated, probably to asymmetric dimethylarginine.,similarity: Contains 1 NTF2 domain.,similarity: Contains 1 RRM (RNA recognition motif) domain.,subunit: Binds to the N-terminal domain of I-kappa-B-alpha.,function: Probable scaffold protein that may be involved in mRNA transport .,PTM: Arg-457 and Arg-468 are dimethylated, probably to asymmetric dimethylarginine.,similarity: Contains 1 NTF2 domain.,similarity: Contains 1 RRM (RNA recognition motif) domain.,subunit: Binds to the N-terminal domain of I-kappa-B-alpha.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from HeLa, HUVEC, and A549 cells, using G3BP2 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HeLa cells using G3BP2 antibody.

产品名称: G3BP2 Rabbit Polyclonal Antibody
产品货号: APRab11213



Western Blot analysis of various cells using G3BP2 Polyclonal Antibody

注意事项 (Note)

For research use only .