

产品名称: Gab 2/4 Rabbit Polyclonal Antibody
产品货号: APRab11225



产品概述 (Summary)

产品名称 (Production Name)	Gab 2/4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GAB2/GAB4 GAB4; GRB2-associated-binding protein 4; GRB2-associated binder 2-like; GAB2-like; GRB2-associated binder 4; GRB2-associated-binding protein 2-like;
别名 (Alternative Names)	Growth factor receptor bound protein 2-associated protein 4; GAB2; KIAA0571; GRB2-associate
基因 ID (Gene ID)	128954.0
蛋白 ID (SwissProt ID)	Q2WGN9/Q9UQC2.The antiserum was produced against synthesized peptide derived from human GAB2. AA range:21-70

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
----------------------	---

产品名称: Gab 2/4 Rabbit Polyclonal Antibody
产品货号: APRab11225




蛋白分子量 (Molecular Weight) 62kDa

研究背景 (Background)

similarity:Belongs to the GAB family.,similarity:Contains 1 PH domain.,similarity:Belongs to the GAB family.,similarity:Contains 1 PH domain.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human tesophagusis, using GAB4 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .