

产品名称: GABAA R $\beta$ 1 Rabbit Polyclonal Antibody  
产品货号: APRab11236



## 产品概述 (Summary)

|                        |                                              |
|------------------------|----------------------------------------------|
| 产品名称 (Production Name) | GABAA R $\beta$ 1 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody                   |
| 宿主 (Host)              | Rabbit                                       |
| 应用 (Application)       | WB,IHC,ICC/IF,ELISA                          |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                              |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                                            |
|------------------------|------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | GABRB1                                                                                                     |
| 别名 (Alternative Names) | GABRB1; Gamma-aminobutyric acid receptor subunit beta-1; GABA(A) receptor subunit beta-1                   |
| 基因 ID (Gene ID)        | 2560.0                                                                                                     |
| 蛋白 ID (SwissProt ID)   | P18505.The antiserum was produced against synthesized peptide derived from human GABA-RB. AA range:401-450 |

## 产品应用 (Application)

|                          |                                                                         |
|--------------------------|-------------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000 |
| 蛋白分子量 (Molecular Weight) | 55kDa                                                                   |

**产品名称: GABAA R $\beta$ 1 Rabbit Polyclonal Antibody**  
**产品货号: AP Rab11236**




## 研究背景 (Background)

The gamma-aminobutyric acid (GABA) A receptor is a multisubunit chloride channel that mediates the fastest inhibitory synaptic transmission in the central nervous system. This gene encodes GABA A receptor, beta 1 subunit. It is mapped to chromosome 4p12 in a cluster comprised of genes encoding alpha 4, alpha 2 and gamma 1 subunits of the GABA A receptor. Alteration of this gene is implicated in the pathogenetics of schizophrenia. [provided by RefSeq, Jul 2008],function:GABA, the major inhibitory neurotransmitter in the vertebrate brain, mediates neuronal inhibition by binding to the GABA/benzodiazepine receptor and opening an integral chloride channel.,online information:Forbidden fruit - Issue 56 of March 2005,similarity:Belongs to the ligand-gated ionic channel (TC 1.A.9) family.,subunit:Binds UBQLN1 (By similarity). Generally pentameric. There are five types of GABA(A) receptor chains: alpha, beta, gamma, delta, and rho.,

## 研究领域 (Research Area)

Neuroactive ligand-receptor interaction;

## 图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using GABA-RB Antibody. The picture on the right is blocked with the synthesized peptide.

## 注意事项 (Note)

For research use only .