

产品名称: GADL1 Rabbit Polyclonal Antibody
产品货号: AP Rab11252



产品概述 (Summary)

产品名称 (Production Name)	GADL1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GADL1
别名 (Alternative Names)	Glutamate decarboxylase-like protein 1 (EC 4.1.1.-)
基因 ID (Gene ID)	339896.0
蛋白 ID (SwissProt ID)	Q6ZQY3.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:200
蛋白分子量 (Molecular Weight)	57kDa

研究背景 (Background)


May catalyze the decarboxylation of L-aspartate, 3-sulfino-L-alanine (cysteine sulfinic acid), and L-cysteate to beta-alanine, hypotaurine and taurine, respectively. Does not exhibit any decarboxylation activity toward glutamate.

产品名称: GADL1 Rabbit Polyclonal Antibody
产品货号: APRab11252



研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of Hela cell lysis, using primary antibody at 1:1000 dilution.

注意事项 (Note)

For research use only .