

产品名称: Gas1 Rabbit Polyclonal Antibody
产品货号: APRab11297



产品概述 (Summary)

产品名称 (Production Name)	Gas1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GAS1
别名 (Alternative Names)	GAS1; Growth arrest-specific protein 1; GAS-1
基因 ID (Gene ID)	2619.0
蛋白 ID (SwissProt ID)	P54826.The antiserum was produced against synthesized peptide derived from human GAS1. AA range:241-290

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	36kDa

研究背景 (Background)

产品名称: Gas1 Rabbit Polyclonal Antibody
产品货号: AP Rab11297




Growth arrest-specific 1 plays a role in growth suppression. GAS1 blocks entry to S phase and prevents cycling of normal and transformed cells. Gas1 is a putative tumor suppressor gene. [provided by RefSeq, Jul 2008],function:Specific growth arrest protein involved in growth suppression. Blocks entry to S phase. Prevents cycling of normal and transformed cells.,


研究领域 (Research Area)

Hedgehog;

图片 (Image Data)



Western blot analysis of lysates from MCF-7 cells, using GAS1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HeLa cells using GAS1 antibody.

注意事项 (Note)

For research use only .