

产品名称: Gasdermin C Rabbit Polyclonal Antibody
产品货号: APRab11304



产品概述 (Summary)

产品名称 (Production Name)	Gasdermin C Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GSDMC
别名 (Alternative Names)	GSDMC; MLZE; Gasdermin-C; Melanoma-derived leucine zipper-containing extranuclear factor
基因 ID (Gene ID)	56169.0
蛋白 ID (SwissProt ID)	Q9BYG8. The antiserum was produced against synthesized peptide derived from human GSDMC. AA range: 109-158

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:100-1:300, ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	55kDa

产品名称: Gasdermin C Rabbit Polyclonal Antibody
产品货号: APRab11304




研究背景 (Background)


caution:Despite its name, it is unclear if this protein contains a functional leucine zipper.,similarity:Belongs to the gasdermin family.,tissue specificity:Expressed primarily in trachea and spleen.,caution:Despite its name, it is unclear if this protein contains a functional leucine zipper.,similarity:Belongs to the gasdermin family.,tissue specificity:Expressed primarily in trachea and spleen.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of GSDMC Antibody. The lane on the right is blocked with the GSDMC peptide.



Western blot analysis of the lysates from HepG2 cells using GSDMC antibody.

产品名称: Gasdermin C Rabbit Polyclonal Antibody
产品货号: APRab11304



Western Blot analysis of various cells using Gasdermin C Polyclonal Antibody diluted at 1: 500

注意事项 (Note)

For research use only .