

产品名称: GCNT6 Rabbit Polyclonal Antibody
产品货号: APRab11362



产品概述 (Summary)

产品名称 (Production Name)	GCNT6 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GCNT6
别名 (Alternative Names)	GCNT6; Beta-1; 3-galactosyl-O-glycosyl-glycoprotein beta-1,6-N-acetylglucosaminyltransferase 6
基因 ID (Gene ID)	644378.0
蛋白 ID (SwissProt ID)	Q5T4J0. The antiserum was produced against synthesized peptide derived from human GCNT6. AA range: 201-250

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	45kDa

产品名称: GCNT6 Rabbit Polyclonal Antibody
产品货号: APRab11362




研究背景 (Background)


function:Glycosyltransferase.,pathway:Protein modification; protein glycosylation.,similarity:Belongs to the glycosyltransferase 14 family.,function:Glycosyltransferase.,pathway:Protein modification; protein glycosylation.,similarity:Belongs to the glycosyltransferase 14 family.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from NIH/3T3 and K562 cells, using GCNT6 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HeLa cells using GCNT6 antibody.

注意事项 (Note)

For research use only .