

产品名称: GCP6 Rabbit Polyclonal Antibody
产品货号: APRab11368



产品概述 (Summary)

产品名称 (Production Name)	GCP6 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TUBGCP6
别名 (Alternative Names)	TUBGCP6; GCP6; KIAA1669; Gamma-tubulin complex component 6; GCP-6
基因 ID (Gene ID)	85378.0
蛋白 ID (SwissProt ID)	Q96RT7.The antiserum was produced against synthesized peptide derived from human TUBGCP6. AA range:741-790

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	200kDa

研究背景 (Background)


产品名称: GCP6 Rabbit Polyclonal Antibody
产品货号: AP Rab11368




The protein encoded by this gene is part of a large multisubunit complex required for microtubule nucleation at the centrosome. [provided by RefSeq, Jul 2008],function:Gamma-tubulin complex is necessary for microtubule nucleation at the centrosome.,sequence caution:Intron retention.,similarity:Belongs to the GCP family.,subunit:Gamma-tubulin complex is composed of gamma-tubulin, GCP2, GCP3, GCP4, GCP5 and GCP6.,

研究领域 (Research Area)

图片 (Image Data)




Immunohistochemistry analysis of paraffin-embedded human brain tissue, using TUBGCP6 Antibody. The picture on the right is blocked with the synthesized peptide.




Western blot analysis of lysates from COLO cells, using TUBGCP6 Antibody. The lane on the right is blocked with the synthesized peptide.


产品名称: GCP6 Rabbit Polyclonal Antibody
产品货号: APRab11368



Western blot analysis of the lysates from HeLa cells using TUBGCP6 antibody.



Western Blot analysis of various cells using GCP6 Polyclonal Antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded Human brain. Antibody was diluted at 1:100 (4°, overnight) . High-pressure and temperature Tris-EDTA, pH 8.0 was used for antigen retrieval. Negative control (right) obtained from antibody was pre-absorbed by immunogen peptide.

注意事项 (Note)

For research use only .