

产品名称: Ghrelin Rabbit Polyclonal Antibody
产品货号: APRab11438



产品概述 (Summary)

产品名称 (Production Name)	Ghrelin Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GHRL
别名 (Alternative Names)	GHRL; MTLRP; Appetite-regulating hormone; Growth hormone secretagogue; Growth hormone-releasing peptide; Motilin-related peptide; Protein M46
基因 ID (Gene ID)	51738.0
蛋白 ID (SwissProt ID)	Q9UBU3.The antiserum was produced against synthesized peptide derived from human Ghrelin. AA range:47-96

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:100-1:300,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	13kDa

产品名称: Ghrelin Rabbit Polyclonal Antibody
产品货号: AP Rab11438




研究背景 (Background)

This gene encodes the ghrelin-obestatin preproprotein that is cleaved to yield two peptides, ghrelin and obestatin. Ghrelin is a powerful appetite stimulant and plays an important role in energy homeostasis. Its secretion is initiated when the stomach is empty, whereupon it binds to the growth hormone secretagogue receptor in the hypothalamus which results in the secretion of growth hormone (somatotropin). Ghrelin is thought to regulate multiple activities, including hunger, reward perception via the mesolimbic pathway, gastric acid secretion, gastrointestinal motility, and pancreatic glucose-stimulated insulin secretion. It was initially proposed that obestatin plays an opposing role to ghrelin by promoting satiety and thus decreasing food intake, but this action is still debated. Recent reports suggest multiple metabolic roles for obestatin, including regulating adipocyte function. Ghrelin is the ligand for growth hormone secretagogue receptor type 1 (GHSR). Induces the release of growth hormone from the pituitary. Has an appetite-stimulating effect, induces adiposity and stimulates gastric acid secretion. Involved in growth regulation. function: Obestatin may be the ligand for GPR39. May have an appetite-reducing effect resulting in decreased food intake. May reduce gastric emptying activity and jejunal motility. mass spectrometry: Ghrelin-27-C10, O-decanoylated form PubMed:12414809, mass spectrometry: Ghrelin-27-C8, O-octanoylated form PubMed:12414809, mass spectrometry: Ghrelin-28-C10, O-decanoylated form PubMed:12414809, mass spectrometry: Ghrelin-28-C10:1, O-decanoylated form PubMed:12414809, mass spectrometry: Ghrelin-28-C8, O-octanoylated form PubMed:12414809, online information: Ghrelin entry, online information: Gut feelings - Issue 66 of January 2006, PTM: Amidation of Leu-98 is essential for obestatin activity, PTM: O-octanoylation or O-decanoylation is essential for ghrelin activity. The O-decanoylated forms Ghrelin-27-C10 and Ghrelin-28-C10 differ in the length of the carbon backbone of the carboxylic acid bound to Ser-26. A small fraction of ghrelin, ghrelin-28-C10:1, may be modified with a singly unsaturated carboxylic acid. similarity: Belongs to the motilin family. tissue specificity: Highest level in stomach. All forms are found in serum as well. Other tissues compensate for the loss of ghrelin synthesis in the stomach following gastrectomy.

研究领域 (Research Area)

图片 (Image Data)




产品名称: Ghrelin Rabbit Polyclonal Antibody
产品货号: AP Rab11438




Immunohistochemistry analysis of Ghrelin antibody in paraffin-embedded human lung carcinoma tissue.



Western blot analysis of lysate from 293 cells, using Ghrelin antibody.




Western Blot analysis of various cells using Ghrelin Polyclonal Antibody



Western Blot analysis of 293 cells using Ghrelin Polyclonal Antibody

产品名称: Ghrelin Rabbit Polyclonal Antibody
产品货号: APRab11438



Immunohistochemical analysis of paraffin-embedded Human lung cancer. Antibody was diluted at 1:100 (4°,overnight) . High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negative contrl (right) obtaned from antibody was pre-absorbed by immunogen peptide.

注意事项 (Note)

For research use only .