

产品名称: GIMAP2 Rabbit Polyclonal Antibody
产品货号: APRab11442



产品概述 (Summary)

产品名称 (Production Name)	GIMAP2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GIMAP2
别名 (Alternative Names)	GIMAP2; IMAP2; GTPase IMAP family member 2; Immunity-associated protein 2; hIMAP2
基因 ID (Gene ID)	26157.0
蛋白 ID (SwissProt ID)	Q9UG22. The antiserum was produced against synthesized peptide derived from human GIMAP2. AA range: 201-250

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	38kDa


产品名称: GIMAP2 Rabbit Polyclonal Antibody
产品货号: APRab11442

研究背景 (Background)


This gene encodes a protein belonging to the GTP-binding superfamily and to the immuno-associated nucleotide (IAN) subfamily of nucleotide-binding proteins. In humans, the IAN subfamily genes are located in a cluster at 7q36.1. [provided by RefSeq, Jul 2008],similarity:Belongs to the IAN GTP-binding protein family.,

研究领域 (Research Area)

图片 (Image Data)




Immunofluorescence analysis of A549 cells, using GIMAP2 Antibody. The picture on the right is blocked with the synthesized



Western blot analysis of lysates from HT-29 cells, using GIMAP2 Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: GIMAP2 Rabbit Polyclonal Antibody
产品货号: APRab11442



Western Blot analysis of various cells using GIMAP2 Polyclonal Antibody

注意事项 (Note)

For research use only .