

产品名称: Glypican-5 Rabbit Polyclonal Antibody
产品货号: AP Rab11523



产品概述 (Summary)

产品名称 (Production Name)	Glypican-5 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPC5
别名 (Alternative Names)	GPC5; Glypican-5
基因 ID (Gene ID)	2262.0
蛋白 ID (SwissProt ID)	P78333.Synthesized peptide derived from the Internal region of human Glypican-5.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	65kDa

研究背景 (Background)


产品名称: Glypican-5 Rabbit Polyclonal Antibody
产品货号: APRab11523



Cell surface heparan sulfate proteoglycans are composed of a membrane-associated protein core substituted with a variable number of heparan sulfate chains. Members of the glypican-related integral membrane proteoglycan family (GRIPS) contain a core protein anchored to the cytoplasmic membrane via a glycosyl phosphatidylinositol linkage. These proteins may play a role in the control of cell division and growth regulation. [provided by RefSeq, Jul 2008],function:Cell surface proteoglycan that bears heparan sulfate.,similarity:Belongs to the glypican family.,tissue specificity:In adult, primarily expressed in the brain. Also detected in fetal brain, lung and liver.,

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of K562 cells using Glypican-5 Polyclonal Antibody

注意事项 (Note)

For research use only .