

产品名称: GM2/GD2 synthase Rabbit Polyclonal Antibody
产品货号: APRab11529



产品概述 (Summary)

产品名称 (Production Name)	GM2/GD2 synthase Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	B4GALNT1
别名 (Alternative Names)	B4GALNT1; GALGT; SIAT2; Beta-1; 4 N-acetylgalactosaminyltransferase 1; (N-acetylneuraminy)-galactosylglucosylceramide; GM2/GD2 synthase; GalNAc-T
基因 ID (Gene ID)	2583.0
蛋白 ID (SwissProt ID)	Q00973.The antiserum was produced against synthesized peptide derived from human B4GALNT1. AA range:281-330

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	60kDa

产品名称: GM2/GD2 synthase Rabbit Polyclonal Antibody
产品货号: APRab11529




研究背景 (Background)

beta-1,4-N-acetyl-galactosaminyltransferase 1(B4GALNT1) Homo sapiens GM2 and GD2 gangliosides are sialic acid-containing glycosphingolipids. GalNAc-T is the enzyme involved in the biosynthesis of G(M2) and G(D2) glycosphingolipids. GalNAc-T catalyzes the transfer of GalNAc into G(M3) and G(D3) by a beta-1,4 linkage, resulting in the synthesis of G(M2) and G(D2), respectively. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2013],catalytic activity:UDP-N-acetyl-D-galactosamine + 1-O-(O-(N-acetyl-alpha-neuraminosyl)-(2->3)-O-beta-D-galactopyranosyl-(1->4)-beta-D-glucopyranosyl)-ceramide = UDP + 1-O-(O-2-(acetylamino)-2-deoxy-beta-D-galactopyranosyl-(1->4)-O-(N-acetyl-alpha-neuraminosyl-(2->3))-O-beta-D-galactopyranosyl-(1->4)-beta-D-glucopyranosyl)-ceramide.,function:Involved in the biosynthesis of gangliosides GM2, GD2 and GA2.,online information:Beta-1,4 N-acetylgalactosaminyltransferase 1,online information:GlycoGene database,pathway:Protein modification; protein glycosylation.,similarity:Belongs to the glycosyltransferase 2 family.,subunit:Homodimer; disulfide-linked.,


研究领域 (Research Area)

Glycosphingolipid biosynthesis;

图片 (Image Data)



Western blot analysis of the lysates from COLO205 cells using B4GALNT1 antibody.



Western Blot analysis of various cells using GM2/GD2 synthase Polyclonal Antibody

产品名称: GM2/GD2 synthase Rabbit Polyclonal Antibody
产品货号: APRab11529



注意事项 (Note)

For research use only .