

产品名称: GOT1B Rabbit Polyclonal Antibody
产品货号: APRab11584



产品概述 (Summary)

产品名称 (Production Name)	GOT1B Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GOLT1B GOT1A CGI-141 HDCMA39P UNQ432/PRO793
别名 (Alternative Names)	
基因 ID (Gene ID)	51026.0
蛋白 ID (SwissProt ID)	Q9Y3E0.Synthesized peptide derived from human protein . at AA range: 70-150

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	15kDa

研究背景 (Background)

function:May be involved in fusion of ER-derived transport vesicles with the Golgi complex.,similarity:Belongs to the GOT1


产品名称: GOT1B Rabbit Polyclonal Antibody
产品货号: APRab11584



family, function: May be involved in fusion of ER-derived transport vesicles with the Golgi complex, similarity: Belongs to the GOT1 family.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night

注意事项 (Note)

For research use only.