

产品名称: GP141 Rabbit Polyclonal Antibody
产品货号: APRab11590



产品概述 (Summary)

产品名称 (Production Name)	GP141 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPR141 PGR13
别名 (Alternative Names)	
基因 ID (Gene ID)	353345.0
蛋白 ID (SwissProt ID)	Q7Z602.Synthesized peptide derived from part region of human protein

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)

GPR141 is a member of the rhodopsin family of G protein-coupled receptors (GPRs) (Fredriksson et al., 2003 [PubMed 14623098]).[supplied by OMIM, Mar 2008],function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 1


产品名称: GP141 Rabbit Polyclonal Antibody
产品货号: APRab11590



family.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from HCT116 cells, primary antibody was diluted at 1:1000, 4° over night

注意事项 (Note)

For research use only .