

产品名称: GPR101 Rabbit Polyclonal Antibody
产品货号: AP Rab11620



产品概述 (Summary)

产品名称 (Production Name)	GPR101 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Mouse, Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPR101
别名 (Alternative Names)	GPR101; Probable G-protein coupled receptor 101
基因 ID (Gene ID)	83550.0
蛋白 ID (SwissProt ID)	Q96P66. The antiserum was produced against synthesized peptide derived from human GPR101. AA range: 451-500

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	60kDa

研究背景 (Background)


产品名称: GPR101 Rabbit Polyclonal Antibody
产品货号: AP Rab11620




The protein encoded by this gene is an orphan G protein-coupled receptor of unknown function. The encoded protein is a member of a family of proteins that contain seven transmembrane domains and transduce extracellular signals through heterotrimeric G proteins. [provided by RefSeq, Sep 2011],function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 1 family.,

研究领域 (Research Area)

图片 (Image Data)




Immunofluorescence analysis of HeLa cells, using GPR101 Antibody. The picture on the right is blocked with the synthesized



Western blot analysis of lysates from COS7 cells, using GPR101 Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: GPR101 Rabbit Polyclonal Antibody
产品货号: AP Rab11620



Western Blot analysis of various cells using GPR101 Polyclonal Antibody diluted at 1: 500

注意事项 (Note)

For research use only .