

产品名称: GPR107 Rabbit Polyclonal Antibody
产品货号: APRab11623



产品概述 (Summary)

产品名称 (Production Name)	GPR107 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPR107
别名 (Alternative Names)	GPR107; KIAA1624; LUSTR1; Protein GPR107; Lung seven transmembrane receptor 1
基因 ID (Gene ID)	57720.0
蛋白 ID (SwissProt ID)	Q5VW38.The antiserum was produced against synthesized peptide derived from human GPR107. AA range:141-190

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	67kDa

产品名称: GPR107 Rabbit Polyclonal Antibody
产品货号: AP Rab11623




研究背景 (Background)

similarity:Belongs to the LU7TM family.,similarity:Belongs to the LU7TM family.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from HUVEC, Jurkat, and COLO205 cells, using GPR107 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using GPR107 Polyclonal Antibody

注意事项 (Note)

For research use only .