

产品名称: GPR108 Rabbit Polyclonal Antibody
产品货号: AP Rab11624



产品概述 (Summary)

产品名称 (Production Name)	GPR108 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPR108
别名 (Alternative Names)	GPR108; LUSTR2; Protein GPR108; Lung seven transmembrane receptor 2
基因 ID (Gene ID)	56927.0
蛋白 ID (SwissProt ID)	Q9NPR9. The antiserum was produced against synthesized peptide derived from human GPR108. AA range: 155-204

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	60kDa

研究背景 (Background)


产品名称: GPR108 Rabbit Polyclonal Antibody
产品货号: APRab11624




similarity:Belongs to the LU7TM family.,similarity:Belongs to the LU7TM family.,

研究领域 (Research Area)


图片 (Image Data)



Immunofluorescence analysis of HUVEC cells, using GPR108 Antibody. The picture on the right is blocked with the synthesized peptide.




Western blot analysis of lysates from A549 cells, using GPR108 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HT-29 cells using GPR108 antibody.

产品名称: GPR108 Rabbit Polyclonal Antibody
产品货号: APRab11624



Western Blot analysis of various cells using GPR108 Polyclonal Antibody

注意事项 (Note)

For research use only .