

产品名称: GPR110 Rabbit Polyclonal Antibody
产品货号: AP Rab11625



产品概述 (Summary)

产品名称 (Production Name)	GPR110 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPR110
别名 (Alternative Names)	GPR110; PGR19; Probable G-protein coupled receptor 110; G-protein coupled receptor KPG_012; G-protein coupled receptor PGR19
基因 ID (Gene ID)	266977.0
蛋白 ID (SwissProt ID)	Q5T601. The antiserum was produced against synthesized peptide derived from human GPR110. AA range: 831-880

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	110kDa

产品名称: GPR110 Rabbit Polyclonal Antibody
产品货号: APRab11625




研究背景 (Background)

function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 2 family. LN-TM7


subfamily.,similarity:Contains 1 GPS domain.,similarity:Contains 1 SEA domain.,function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 2 family. LN-TM7 subfamily.,similarity:Contains 1 GPS domain.,similarity:Contains 1 SEA domain.,

研究领域 (Research Area)

图片 (Image Data)




Immunofluorescence analysis of A549 cells, using GPR110 Antibody. The picture on the right is blocked with the synthesized




Western blot analysis of lysates from Jurkat, HUVEC, and COLO cells, using GPR110 Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: GPR110 Rabbit Polyclonal Antibody
产品货号: AP Rab11625



Western Blot analysis of various cells using GPR110 Polyclonal Antibody



Western Blot analysis of 3T3 cells using GPR110 Polyclonal Antibody

注意事项 (Note)

For research use only .