

产品名称: GPR113 Rabbit Polyclonal Antibody
产品货号: AP Rab11626



产品概述 (Summary)

产品名称 (Production Name)	GPR113 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPR113
别名 (Alternative Names)	GPR113; PGR23; Probable G-protein coupled receptor 113; G-protein coupled receptor PGR23
基因 ID (Gene ID)	165082.0
蛋白 ID (SwissProt ID)	Q8IZF5.The antiserum was produced against synthesized peptide derived from human GPR113. AA range:21-70

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	116kDa

产品名称: GPR113 Rabbit Polyclonal Antibody
产品货号: AP Rab11626




研究背景 (Background)

function: Orphan receptor., similarity: Belongs to the G-protein coupled receptor 2 family. LN-TM7


subfamily., similarity: Contains 1 GPS domain., function: Orphan receptor., similarity: Belongs to the G-protein coupled receptor 2 family. LN-TM7 subfamily., similarity: Contains 1 GPS domain.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from LOVO cells, using GPR113 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of HeLa cells using GPR113 Polyclonal Antibody diluted at 1: 2000

注意事项 (Note)

For research use only .