

产品名称: GPR114 Rabbit Polyclonal Antibody
产品货号: APRab11627



产品概述 (Summary)

产品名称 (Production Name)	GPR114 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPR114
别名 (Alternative Names)	GPR114; PGR27; Probable G-protein coupled receptor 114; G-protein coupled receptor PGR27
基因 ID (Gene ID)	221188.0
蛋白 ID (SwissProt ID)	Q8IZF4. The antiserum was produced against synthesized peptide derived from human GPR114. AA range: 91-140

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	59kDa

产品名称: GPR114 Rabbit Polyclonal Antibody
产品货号: AP Rab11627




研究背景 (Background)


This gene encodes a member of the adhesion family of G-protein coupled receptors. Members of this family are characterized by long N-termini and multiple functional domains. They may play a role in the immune system as well as in the central nervous system. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2016],function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 2 family. LN-TM7 subfamily.,similarity:Contains 1 GPS domain.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from COLO205 cells, using GPR114 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from 293 cells using GPR114 antibody.

产品名称: GPR114 Rabbit Polyclonal Antibody
产品货号: APRab11627



Western Blot analysis of COLO205 cells using GPR114 Polyclonal Antibody

注意事项 (Note)

For research use only .