

**产品名称: GPR115 Rabbit Polyclonal Antibody**  
**产品货号: APRab11628**



## 产品概述 (Summary)

|                        |                                   |
|------------------------|-----------------------------------|
| 产品名称 (Production Name) | GPR115 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody        |
| 宿主 (Host)              | Rabbit                            |
| 应用 (Application)       | WB, ICC/IF, ELISA                 |
| 种属反应性 (Reactivity)     | Human, Mouse                      |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                                             |
|------------------------|-------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | GPR115                                                                                                      |
| 别名 (Alternative Names) | GPR115; PGR18; Probable G-protein coupled receptor 115; G-protein coupled receptor PGR18                    |
| 基因 ID (Gene ID)        | 221393.0                                                                                                    |
| 蛋白 ID (SwissProt ID)   | Q8IZF3. The antiserum was produced against synthesized peptide derived from human GPR115. AA range: 641-690 |

## 产品应用 (Application)

|                          |                                                            |
|--------------------------|------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:5000-1:20000 |
| 蛋白分子量 (Molecular Weight) | 82kDa                                                      |

产品名称: GPR115 Rabbit Polyclonal Antibody  
产品货号: AP Rab11628




## 研究背景 (Background)

function: Orphan receptor., similarity: Belongs to the G-protein coupled receptor 2 family. LN-TM7


subfamily., similarity: Contains 1 GPS domain., function: Orphan receptor., similarity: Belongs to the G-protein coupled receptor 2 family. LN-TM7 subfamily., similarity: Contains 1 GPS domain.,

## 研究领域 (Research Area)

## 图片 (Image Data)




Immunofluorescence analysis of LOVO cells, using GPR115 Antibody. The picture on the right is blocked with the synthesized



Western blot analysis of lysates from RAW264.7 cells, using GPR115 Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: GPR115 Rabbit Polyclonal Antibody  
产品货号: AP Rab11628



Western Blot analysis of RAW264.7 cells using GPR115 Polyclonal Antibody

### 注意事项 (Note)

For research use only .