

产品名称: GPR144 Rabbit Polyclonal Antibody
产品货号: AP Rab11643



产品概述 (Summary)

产品名称 (Production Name)	GPR144 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPR144
别名 (Alternative Names)	GPR144; PGR24; Probable G-protein coupled receptor 144; G-protein coupled receptor PGR24
基因 ID (Gene ID)	347088.0
蛋白 ID (SwissProt ID)	Q7Z7M1. The antiserum was produced against synthesized peptide derived from human GPR144. AA range: 91-140

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	104kDa

产品名称: GPR144 Rabbit Polyclonal Antibody
产品货号: APRab11643




研究背景 (Background)

function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 2 family. LN-TM7
subfamily.,similarity:Contains 1 GPS domain.,similarity:Contains 1 pentaxin domain.,function:Orphan
receptor.,similarity:Belongs to the G-protein coupled receptor 2 family. LN-TM7 subfamily.,similarity:Contains 1 GPS
domain.,similarity:Contains 1 pentaxin domain.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from K562 cells, using GPR144 Antibody. The lane on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .