

产品名称: GPR15 Rabbit Polyclonal Antibody
产品货号: AP Rab11646



产品概述 (Summary)

产品名称 (Production Name)	GPR15 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPR15
别名 (Alternative Names)	GPR15; G-protein coupled receptor 15; Brother of Bonzo; BoB
基因 ID (Gene ID)	2838.0
蛋白 ID (SwissProt ID)	P49685. The antiserum was produced against synthesized peptide derived from human GPR15. AA range: 201-250

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	40kDa

研究背景 (Background)


产品名称: GPR15 Rabbit Polyclonal Antibody
产品货号: AP Rab11646




This gene encodes a G protein-coupled receptor that acts as a chemokine receptor for human immunodeficiency virus type 1 and 2. The encoded protein localizes to the cell membrane. [provided by RefSeq, Nov 2012],function:Probable chemokine receptor. Alternative coreceptor with CD4 for HIV-1 infection.,similarity:Belongs to the G-protein coupled receptor 1 family.,

研究领域 (Research Area)

图片 (Image Data)




Immunofluorescence analysis of LOVO cells, using GPR15 Antibody. The picture on the right is blocked with the synthesized




Western blot analysis of lysates from HepG2 cells, using GPR15 Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: GPR15 Rabbit Polyclonal Antibody
产品货号: AP Rab11646



Western blot analysis of the lysates from HeLa cells using GPR15 antibody.



Western Blot analysis of various cells using GPR15 Polyclonal Antibody

注意事项 (Note)

For research use only .