

产品名称: GPR153 Rabbit Polyclonal Antibody
产品货号: AP Rab11650



产品概述 (Summary)

产品名称 (Production Name)	GPR153 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPR153
别名 (Alternative Names)	GPR153; PGR1; Probable G-protein coupled receptor 153; G-protein coupled receptor PGR1
基因 ID (Gene ID)	387509.0
蛋白 ID (SwissProt ID)	Q6NV75.The antiserum was produced against synthesized peptide derived from human GPR153. AA range:195-244

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	65kDa

产品名称: GPR153 Rabbit Polyclonal Antibody
产品货号: AP Rab11650




研究背景 (Background)


function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 1 family.,function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 1 family.,

研究领域 (Research Area)

图片 (Image Data)




Immunofluorescence analysis of LOVO cells, using GPR153 Antibody. The picture on the right is blocked with the synthesized



Western blot analysis of lysates from K562 cells, using GPR153 Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: GPR153 Rabbit Polyclonal Antibody
产品货号: APRab11650



Western blot analysis of the lysates from HepG2 cells using GPR153 antibody.

注意事项 (Note)

For research use only .