

产品名称: GPR171 Rabbit Polyclonal Antibody
产品货号: APRab11656



产品概述 (Summary)

产品名称 (Production Name)	GPR171 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPR171
别名 (Alternative Names)	GPR171; H963; Probable G-protein coupled receptor 171; G-protein coupled receptor H963
基因 ID (Gene ID)	29909.0
蛋白 ID (SwissProt ID)	O14626. The antiserum was produced against synthesized peptide derived from human GPR171. AA range: 181-230

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	36kDa

产品名称: GPR171 Rabbit Polyclonal Antibody
产品货号: AP Rab11656




研究背景 (Background)


function: Orphan receptor., similarity: Belongs to the G-protein coupled receptor 1 family., function: Orphan receptor., similarity: Belongs to the G-protein coupled receptor 1 family.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from HT-29 cells, using GPR171 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using GPR171 Polyclonal Antibody

注意事项 (Note)

For research use only .