

产品名称: GPR175 Rabbit Polyclonal Antibody
产品货号: AP Rab11661



产品概述 (Summary)

产品名称 (Production Name)	GPR175 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TPRA1
别名 (Alternative Names)	TPRA1; GPR175; PP6566; Transmembrane protein adipocyte-associated 1; Integral membrane protein GPR175
基因 ID (Gene ID)	131601.0
蛋白 ID (SwissProt ID)	Q86W33. The antiserum was produced against synthesized peptide derived from human GPR175. AA range: 198-247

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	41kDa

产品名称: GPR175 Rabbit Polyclonal Antibody
产品货号: AP Rab11661




研究背景 (Background)


TPRA1 is a gene on chromosome 3q21.3 that encodes transmembrane protein adipocyte-associated 1 expressed in most tissues, especially in heart, placenta and kidney, and plays role in lipid metabolism and in ageing. Gene Ontology (GO) annotations related to this gene include G-protein coupled receptor activity. aging,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from 293 cells, using GPR175 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using GPR175 Polyclonal Antibody diluted at 1: 2000

注意事项 (Note)

For research use only .