

产品名称: GPR180 Rabbit Polyclonal Antibody
产品货号: AP Rab11664



产品概述 (Summary)

产品名称 (Production Name)	GPR180 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPR180
别名 (Alternative Names)	GPR180; ITR; Integral membrane protein GPR180; Intimal thickness-related receptor
基因 ID (Gene ID)	160897.0
蛋白 ID (SwissProt ID)	Q86V85.The antiserum was produced against synthesized peptide derived from human GPR180. AA range:241-290

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	50kDa

产品名称: GPR180 Rabbit Polyclonal Antibody
产品货号: AP Rab11664




研究背景 (Background)


This gene encodes a protein that is a member of the G protein-coupled receptor superfamily. This protein is produced predominantly in vascular smooth muscle cells and may play an important role in the regulation of vascular remodeling. [provided by RefSeq, Jul 2008],

研究领域 (Research Area)

图片 (Image Data)




Immunofluorescence analysis of A549 cells, using GPR180 Antibody. The picture on the right is blocked with the synthesized peptide.




Immunohistochemistry analysis of paraffin-embedded human brain tissue, using GPR180 Antibody. The picture on the right is blocked with the synthesized peptide.

产品名称: GPR180 Rabbit Polyclonal Antibody
产品货号: AP Rab11664



Western blot analysis of the lysates from Jurkat cells using GPR180 antibody.



Western Blot analysis of LOVO cells using GPR180 Polyclonal Antibody

注意事项 (Note)

For research use only .