

产品名称: GPR19 Rabbit Polyclonal Antibody
产品货号: APRab11665



产品概述 (Summary)

产品名称 (Production Name)	GPR19 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPR19
别名 (Alternative Names)	GPR19; Probable G-protein coupled receptor 19; GPR-NGA
基因 ID (Gene ID)	2842.0
蛋白 ID (SwissProt ID)	Q15760.The antiserum was produced against synthesized peptide derived from human GPR19. AA range:361-410

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	47kDa

研究背景 (Background)


产品名称: GPR19 Rabbit Polyclonal Antibody
产品货号: AP Rab11665




function: Orphan receptor., similarity: Belongs to the G-protein coupled receptor 1 family., tissue specificity: Abundant expression in the brain., function: Orphan receptor., similarity: Belongs to the G-protein coupled receptor 1 family., tissue specificity: Abundant expression in the brain.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from MCF-7, HUVEC, and HepG2 cells, using GPR19 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using GPR19 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA) .

注意事项 (Note)

For research use only .