

产品名称: GPR32 Rabbit Polyclonal Antibody
产品货号: APRab11674



产品概述 (Summary)

产品名称 (Production Name)	GPR32 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPR32
别名 (Alternative Names)	GPR32; Probable G-protein coupled receptor 32
基因 ID (Gene ID)	2854.0
蛋白 ID (SwissProt ID)	O75388. The antiserum was produced against synthesized peptide derived from human GPR32. AA range: 151-200

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	33kDa

研究背景 (Background)


产品名称: GPR32 Rabbit Polyclonal Antibody
产品货号: AP Rab11674




This gene is intronless and encodes a member of the G-protein coupled receptor 1 family. The encoded protein binds to resolvin D1 and lipoxin A4 and has been linked to pulmonary inflammation. A related pseudogene has been identified on chromosome 19. [provided by RefSeq, Nov 2012],function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 1 family.,

研究领域 (Research Area)

图片 (Image Data)



Immunofluorescence analysis of A549 cells, using GPR32 Antibody. The picture on the right is blocked with the synthesized



Western blot analysis of lysates from A549 cells, using GPR32 Antibody. The lane on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .