

产品名称: GPR34 Rabbit Polyclonal Antibody
产品货号: APRab11676



产品概述 (Summary)

产品名称 (Production Name)	GPR34 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Mouse, Rat, Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPR34
别名 (Alternative Names)	GPR34; Probable G-protein coupled receptor 34
基因 ID (Gene ID)	2857.0
蛋白 ID (SwissProt ID)	Q9UPC5. The antiserum was produced against synthesized peptide derived from human GPR34. AA range: 181-230

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	44kDa

研究背景 (Background)


产品名称: GPR34 Rabbit Polyclonal Antibody
产品货号: AP Rab11676




G protein-coupled receptors (GPCRs), such as GPR34, are integral membrane proteins containing 7 putative transmembrane domains (TMs). These proteins mediate signals to the interior of the cell via activation of heterotrimeric G proteins that in turn activate various effector proteins, ultimately resulting in a physiologic response.[supplied by OMIM, Apr 2006],function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 1 family.,tissue specificity:Broadly expressed.,

研究领域 (Research Area)

图片 (Image Data)




Immunofluorescence analysis of LOVO cells, using GPR34 Antibody. The picture on the right is blocked with the synthesized



Western blot analysis of lysates from HT-29, K562, and COS7 cells, using GPR34 Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: GPR34 Rabbit Polyclonal Antibody
产品货号: AP Rab11676



Western Blot analysis of various cells using GPR34 Polyclonal Antibody

注意事项 (Note)

For research use only .