

产品名称: GPR38 Rabbit Polyclonal Antibody
产品货号: APRab11680



产品概述 (Summary)

产品名称 (Production Name)	GPR38 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MLNR
别名 (Alternative Names)	MLNR; GPR38; MTLR; MTLR1; Motilin receptor; G-protein coupled receptor 38
基因 ID (Gene ID)	2862.0
蛋白 ID (SwissProt ID)	O43193.Synthesized peptide derived from the C-terminal region of human GPR38.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	45kDa

研究背景 (Background)

产品名称: GPR38 Rabbit Polyclonal Antibody
产品货号: AP Rab11680




Motilin is a 22 amino acid peptide hormone expressed throughout the gastrointestinal (GI) tract. The protein encoded by this gene is a motilin receptor which is a member of the G-protein coupled receptor 1 family. This member is a multi-pass transmembrane protein, and is an important therapeutic target for the treatment of hypomotility disorders. [provided by RefSeq, Aug 2011],function:Receptor for motilin.,similarity:Belongs to the G-protein coupled receptor 1 family.,tissue specificity:Expressed only in thyroid, stomach, and bone marrow.,

研究领域 (Research Area)

Neuroactive ligand-receptor interaction;

图片 (Image Data)



Western Blot analysis of various cells using GPR38 Polyclonal Antibody

注意事项 (Note)

For research use only .