

产品名称: GPR42 Rabbit Polyclonal Antibody
产品货号: APRab11685



产品概述 (Summary)

产品名称 (Production Name)	GPR42 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPR42
别名 (Alternative Names)	GPR42; GPR42P; G-protein coupled receptor 42
基因 ID (Gene ID)	2866.0
蛋白 ID (SwissProt ID)	O15529.The antiserum was produced against synthesized peptide derived from human GPR42. AA range:12-61

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC/IF 1:200-1:1000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)


产品名称: GPR42 Rabbit Polyclonal Antibody
产品货号: AP Rab11685



GPR42, also named as GPR42P, FFAR1L and GPR41L, belongs to the G-protein coupled receptor 1 family. GPR42 appears to have arisen by tandem duplication of FFAR3 and have lost the ability to activate G-protein family in response to carboxylate ligands primarily due to an amino acid change at position 174. It is an orphan receptor.

研究领域 (Research Area)

图片 (Image Data)



Immunofluorescence analysis of MCF-7 cells, using GPR42 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .