

产品名称: GPR45 Rabbit Polyclonal Antibody
产品货号: AP Rab11687



产品概述 (Summary)

产品名称 (Production Name)	GPR45 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPR45
别名 (Alternative Names)	
基因 ID (Gene ID)	11250.0
蛋白 ID (SwissProt ID)	Q9Y5Y3.Synthesized peptide derived from human protein . at AA range: 110-190

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	40kDa

研究背景 (Background)

This intronless gene encodes a member of the G protein-coupled receptor (GPCR) family. Members of this protein family


产品名称: GPR45 Rabbit Polyclonal Antibody
产品货号: APRab11687



contain seven putative transmembrane domains and may mediate signaling processes to the interior of the cell via activation of heterotrimeric G proteins. This protein may function in the central nervous system. [provided by RefSeq, Jul 2008],function:Orphan receptor. May play a role in brain function.,similarity:Belongs to the G-protein coupled receptor 1 family.,tissue specificity:Expressed in brain; detected in the basal forebrain, frontal cortex, and caudate, but not in thalamus, hippocampus, or putamen.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from Jarkat cells, primary antibody was diluted at 1:1000, 4°over night

注意事项 (Note)

For research use only .