

产品名称: GPR52 Rabbit Polyclonal Antibody
产品货号: APRab11690



产品概述 (Summary)

产品名称 (Production Name)	GPR52 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPR52
别名 (Alternative Names)	GPR52; Probable G-protein coupled receptor 52
基因 ID (Gene ID)	9293.0
蛋白 ID (SwissProt ID)	Q9Y2T5.The antiserum was produced against synthesized peptide derived from human GPR52. AA range:311-360

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	41kDa

研究背景 (Background)


产品名称: GPR52 Rabbit Polyclonal Antibody
产品货号: APRab11690




Members of the G protein-coupled receptor (GPR) family play important roles in signal transduction from the external environment to the inside of the cell.[supplied by OMIM, Jul 2002],function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 1 family.,

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of lysates from COLO, HUVEC, HeLa, and HepG2 cells, using GPR52 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using GPR52 Polyclonal Antibody diluted at 1: 2000



产品名称: GPR52 Rabbit Polyclonal Antibody
产品货号: APRab11690



Western blot analysis of 293T 3T3 lysis using GPR52 antibody. Antibody was diluted at 1:2000

注意事项 (Note)

For research use only .