

产品名称: GPR62 Rabbit Polyclonal Antibody
产品货号: APRab11695



产品概述 (Summary)

产品名称 (Production Name)	GPR62 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPR62
别名 (Alternative Names)	GPR62; Probable G-protein coupled receptor 62; G-protein coupled receptor GPCR8; hGPCR8; G-protein coupled receptor KPG_005
基因 ID (Gene ID)	118442.0
蛋白 ID (SwissProt ID)	Q9BZJ7.The antiserum was produced against synthesized peptide derived from human GPR62. AA range:181-230

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	38kDa

产品名称: GPR62 Rabbit Polyclonal Antibody
产品货号: AP Rab11695




研究背景 (Background)


function: Orphan receptor., similarity: Belongs to the G-protein coupled receptor 1 family., tissue specificity: Expressed in brain; detected in the basal forebrain, frontal cortex, caudate, putamen, thalamus and hippocampus., function: Orphan receptor., similarity: Belongs to the G-protein coupled receptor 1 family., tissue specificity: Expressed in brain; detected in the basal forebrain, frontal cortex, caudate, putamen, thalamus and hippocampus.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from COLO cells, using GPR62 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from Jurkat cells using GPR62 antibody.

产品名称: GPR62 Rabbit Polyclonal Antibody
产品货号: APRab11695



Western Blot analysis of various cells using GPR62 Polyclonal Antibody

注意事项 (Note)

For research use only .