

产品名称: GPR63 Rabbit Polyclonal Antibody
产品货号: APRab11696



产品概述 (Summary)

产品名称 (Production Name)	GPR63 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPR63
别名 (Alternative Names)	GPR63; PSP24B; Probable G-protein coupled receptor 63; PSP24-2; PSP24-beta
基因 ID (Gene ID)	81491.0
蛋白 ID (SwissProt ID)	Q9BZJ6.The antiserum was produced against synthesized peptide derived from human GPR63. AA range:370-419

产品应用 (Application)

稀释比 (Dilution Ratio)	ICC/IF 1:200-1:1000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

产品名称: GPR63 Rabbit Polyclonal Antibody
产品货号: AP Rab11696




研究背景 (Background)

This gene encodes a G protein-coupled receptor. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Dec 2011], function: Orphan receptor. May play a role in brain function., similarity: Belongs to the G-protein coupled receptor 1 family., tissue specificity: Expressed in brain; detected in the frontal cortex, with lower levels in the thalamus, caudate, hypothalamus and midbrain.,

研究领域 (Research Area)

图片 (Image Data)



Immunofluorescence analysis of MCF7 cells, using GPR63 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .