

产品名称: GPR85 Rabbit Polyclonal Antibody
产品货号: APRab11704



产品概述 (Summary)

产品名称 (Production Name)	GPR85 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPR85
别名 (Alternative Names)	GPR85; SREB2; Probable G-protein coupled receptor 85; Super conserved receptor expressed in brain 2
基因 ID (Gene ID)	54329.0
蛋白 ID (SwissProt ID)	P60893. The antiserum was produced against synthesized peptide derived from human GPR85. AA range: 181-230

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:200-1:1000, ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	

产品名称: GPR85 Rabbit Polyclonal Antibody
产品货号: AP Rab11704




研究背景 (Background)

Members of the G protein-coupled receptor (GPCR) family, such as GPR85, have a similar structure characterized by 7 transmembrane domains. Activation of GPCRs by extracellular stimuli, such as neurotransmitters, hormones, or light, induces an intracellular signaling cascade mediated by heterotrimeric GTP-binding proteins, or G proteins (Matsumoto et al., 2000 [PubMed 10833454]).[supplied by OMIM, Aug 2008],function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 1 family.,tissue specificity:Highly expressed in brain and testis. Lower levels in small intestine, placenta and spleen. In brain regions, detected in all regions tested, but somewhat lower levels in the corpus callosum, medulla and spinal cord.,

研究领域 (Research Area)

图片 (Image Data)



Immunofluorescence analysis of LOVO cells, using GPR85 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .