

产品名称: GPR87 Rabbit Polyclonal Antibody
产品货号: APRab11705



产品概述 (Summary)

产品名称 (Production Name)	GPR87 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPR87
别名 (Alternative Names)	GPR87; GPR95; FKSG78; G-protein coupled receptor 87; G-protein coupled receptor 95
基因 ID (Gene ID)	53836.0
蛋白 ID (SwissProt ID)	Q9BY21. The antiserum was produced against synthesized peptide derived from human GPR87. AA range: 221-270

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	40kDa

产品名称: GPR87 Rabbit Polyclonal Antibody
产品货号: AP Rab11705




研究背景 (Background)


This gene encodes a G protein-coupled receptor and is located in a cluster of G protein-coupled receptor genes on chromosome 3. The encoded protein has been shown to be overexpressed in lung squamous cell carcinoma (PMID:18057535) and regulated by p53 (PMID:19602589). [provided by RefSeq, Nov 2011], function: Orphan receptor, similarity: Belongs to the G-protein coupled receptor 1 family, tissue specificity: Expressed in placenta and prostate. Weaker expression in thymus. Not expressed in thalamus, hippocampus, pons or cerebellum.

研究领域 (Research Area)

图片 (Image Data)



Immunofluorescence analysis of HUVEC cells, using GPR87 Antibody. The picture on the right is blocked with the synthesized



Western blot analysis of lysates from HUVEC and HT-29 cells, using GPR87 Antibody. The lane on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .