

产品名称: GPR87/95 Rabbit Polyclonal Antibody
产品货号: APRab11706



产品概述 (Summary)

产品名称 (Production Name)	GPR87/95 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPR87
别名 (Alternative Names)	GPR87; GPR95; FKSG78; G-protein coupled receptor 87; G-protein coupled receptor 95
基因 ID (Gene ID)	53836.0
蛋白 ID (SwissProt ID)	Q9BY21.The antiserum was produced against synthesized peptide derived from the N-terminal region of human GPR87. AA range:1-50

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	42kDa

产品名称: GPR87/95 Rabbit Polyclonal Antibody
产品货号: AP Rab11706




研究背景 (Background)

This gene encodes a G protein-coupled receptor and is located in a cluster of G protein-coupled receptor genes on chromosome 3. The encoded protein has been shown to be overexpressed in lung squamous cell carcinoma (PMID:18057535) and regulated by p53 (PMID:19602589). [provided by RefSeq, Nov 2011], function: Orphan receptor, similarity: Belongs to the G-protein coupled receptor 1 family, tissue specificity: Expressed in placenta and prostate. Weaker expression in thymus. Not expressed in thalamus, hippocampus, pons or cerebellum.

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human-lung, antibody was diluted at 1:100



Western Blot analysis of K562 cells using GPR87/95 Polyclonal Antibody. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only.