

产品名称: GPR92 Rabbit Polyclonal Antibody
产品货号: APRab11709



产品概述 (Summary)

产品名称 (Production Name)	GPR92 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LPAR5
别名 (Alternative Names)	LPAR5; GPR92; GPR93; Lysophosphatidic acid receptor 5; LPA receptor 5; LPA-5; G-protein coupled receptor 92; G-protein coupled receptor 93
基因 ID (Gene ID)	57121.0
蛋白 ID (SwissProt ID)	Q9H1C0.The antiserum was produced against synthesized peptide derived from human GPR92. AA range:241-290

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	40kDa

产品名称: GPR92 Rabbit Polyclonal Antibody
产品货号: AP Rab11709




研究背景 (Background)

lysophosphatidic acid receptor 5(LPAR5) Homo sapiens This gene encodes a member of the rhodopsin class of G protein-coupled transmembrane receptors. This protein transmits extracellular signals from lysophosphatidic acid to cells through heterotrimeric G proteins and mediates numerous cellular processes. Many G protein receptors serve as targets for pharmaceutical drugs. Transcript variants of this gene have been described.[provided by RefSeq, Dec 2008],function:Receptor for lysophosphatidic acid (LPA), a mediator of diverse cellular activities.,similarity:Belongs to the G-protein coupled receptor 1 family.,tissue specificity:Not expressed in frontal cortex, basal forebrain, caudate putamen, thalamus, or hippocampus.,


研究领域 (Research Area)

PI3K/Akt

图片 (Image Data)



Western blot analysis of lysates from HT-29 and K562 cells, using GPR92 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using GPR92 Polyclonal Antibody diluted at 1: 500

注意事项 (Note)

For research use only .

产品名称: GPR92 Rabbit Polyclonal Antibody
产品货号: APRab11709

