

产品名称: Granzyme B/H Rabbit Polyclonal Antibody
产品货号: AP Rab11741



产品概述 (Summary)

产品名称 (Production Name)	Granzyme B/H Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GZMB/GZMH GZMB; CGL1; CSPB; CTLA1; GRB; Granzyme B; C11; CTLA-1; Cathepsin G-like 1;
别名 (Alternative Names)	CTSG1; Cytotoxic T-lymphocyte proteinase 2; Lymphocyte protease; Fragmentin-2; Granzyme-2; Human lymphocyte protein; HLP; SECT; T-cell serine protease 1-3E; GZMH; C
基因 ID (Gene ID)	2999.0
蛋白 ID (SwissProt ID)	P10144/P20718.The antiserum was produced against synthesized peptide derived from human Granzyme B. AA range:10-59

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:10000
----------------------	--

产品名称: Granzyme B/H Rabbit Polyclonal Antibody
产品货号: APRab11741



蛋白分子量 (Molecular Weight) 38kDa


研究背景 (Background)

This gene encodes a member of the granzyme subfamily of proteins, part of the peptidase S1 family of serine proteases. The encoded preproprotein is secreted by natural killer (NK) cells and cytotoxic T lymphocytes (CTLs) and proteolytically processed to generate the active protease, which induces target cell apoptosis. This protein also processes cytokines and degrades extracellular matrix proteins, and these roles are implicated in chronic inflammation and wound healing. Expression of this gene may be elevated in human patients with cardiac fibrosis. [provided by RefSeq, Sep 2016], catalytic activity: Preferential cleavage: -Asp-|-Xaa- >> -Asn-|-Xaa- > -Met-|-Xaa-, -Ser-|-Xaa-, function: This enzyme is necessary for target cell lysis in cell-mediated immune responses. It cleaves after Asp. Seems to be linked to an activation cascade of caspases (aspartate-specific cysteine proteases) responsible for apoptosis execution. Cleaves caspase-3, -7, -9 and 10 to give rise to active enzymes mediating apoptosis, induction: By staphylococcal enterotoxin A (SEA) in peripheral blood leukocytes, similarity: Belongs to the peptidase S1 family, similarity: Belongs to the peptidase S1 family. Granzyme subfamily, similarity: Contains 1 peptidase S1 domain, subcellular location: Cytoplasmic granules of cytolytic T-lymphocytes and natural killer cells,

研究领域 (Research Area)


Natural killer cell mediated cytotoxicity; Type I diabetes mellitus; Autoimmune thyroid disease; Allograft rejection; Graft-versus-host disease;

图片 (Image Data)




Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using Granzyme B Antibody. The picture on the right is blocked with the synthesized peptide.

产品名称: Granzyme B/H Rabbit Polyclonal Antibody
产品货号: APRab11741



Western blot analysis of lysates from NIH/3T3 cells, using Granzyme B Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of 3T3 cells using Granzyme B/H Polyclonal Antibody diluted at 1: 500

注意事项 (Note)

For research use only .