

产品名称: Granzyme M Rabbit Polyclonal Antibody
产品货号: AP Rab11744



产品概述 (Summary)

产品名称 (Production Name)	Granzyme M Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GZMM MET1
别名 (Alternative Names)	Granzyme M (EC 3.4.21.-; Met-1 serine protease; Hu-Met-1; Met-ase; Natural killer cell granular protease)
基因 ID (Gene ID)	3004.0
蛋白 ID (SwissProt ID)	P51124. Synthetic peptide from human protein at AA range: 207-256

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:200, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)


产品名称: Granzyme M Rabbit Polyclonal Antibody
产品货号: APRab11744




Human natural killer (NK) cells and activated lymphocytes express and store a distinct subset of neutral serine proteases together with proteoglycans and other immune effector molecules in large cytoplasmic granules. These serine proteases are collectively termed granzymes and include 4 distinct gene products: granzyme A, granzyme B, granzyme H, and the protein encoded by this gene, granzyme M. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2012],function:Cleaves peptide substrates after methionine, leucine, and norleucine.,similarity:Belongs to the peptidase S1 family. Granzyme subfamily.,similarity:Contains 1 peptidase S1 domain.,subcellular location:Granules of large granular lymphocytes.,

研究领域 (Research Area)

图片 (Image Data)




Immunohistochemical analysis of paraffin-embedded Human-tonsil, antibody was diluted at 1:100



Immunohistochemical analysis of paraffin-embedded Human-lymph, antibody was diluted at 1:100

产品名称: Granzyme M Rabbit Polyclonal Antibody
产品货号: APRab11744



Immunohistochemical analysis of paraffin-embedded Human-lymph, antibody was diluted at 1:100

注意事项 (Note)

For research use only .