

产品名称: GRIF-1 Rabbit Polyclonal Antibody
产品货号: APRab11760



产品概述 (Summary)

产品名称 (Production Name)	GRIF-1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TRAK2
别名 (Alternative Names)	TRAK2; ALS2CR3; KIAA0549; Trafficking kinesin-binding protein 2; Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 3 protein
基因 ID (Gene ID)	66008.0
蛋白 ID (SwissProt ID)	O60296. The antiserum was produced against synthesized peptide derived from human TRAK2. AA range: 430-480

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	

产品名称: GRIF-1 Rabbit Polyclonal Antibody
产品货号: APRab11760




研究背景 (Background)

PTM:O-glycosylated.,similarity:Contains 1 HAP1 N-terminal domain.,subunit:Interacts with GABA-A receptor and O-GlcNAc transferase (By similarity). Interacts with RHOT1/Miro-1 and RHOT2/Miro-2.,tissue specificity:Widely expressed, with highest expression in heart.,PTM:O-glycosylated.,similarity:Contains 1 HAP1 N-terminal domain.,subunit:Interacts with GABA-A receptor and O-GlcNAc transferase (By similarity). Interacts with RHOT1/Miro-1 and RHOT2/Miro-2.,tissue specificity:Widely expressed, with highest expression in heart.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200 (4° overnight) . 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200 (room temperature, 30min) .

注意事项 (Note)

For research use only .