

产品名称: GRIN2 Rabbit Polyclonal Antibody
产品货号: APRab11766



产品概述 (Summary)

产品名称 (Production Name)	GRIN2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GPRIN2
别名 (Alternative Names)	GPRIN2; KIAA0514; G protein-regulated inducer of neurite outgrowth 2; GRIN2
基因 ID (Gene ID)	9721.0
蛋白 ID (SwissProt ID)	O60269.The antiserum was produced against synthesized peptide derived from human GPRIN2. AA range:11-60

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	47kDa

产品名称: GRIN2 Rabbit Polyclonal Antibody
产品货号: APRab11766




研究背景 (Background)


function: May be involved in neurite outgrowth., subunit: Interacts with activated forms of GNAO1 and GNAZ., tissue specificity: Expressed specifically in the cerebellum., function: May be involved in neurite outgrowth., subunit: Interacts with activated forms of GNAO1 and GNAZ., tissue specificity: Expressed specifically in the cerebellum.,

研究领域 (Research Area)

图片 (Image Data)




Western Blot analysis of various cells using GRIN2 Polyclonal Antibody



Immunohistochemical analysis of paraffin-embedded Human brain. Antibody was diluted at 1:100 (4°, overnight) . High-pressure and temperature Tris-EDTA, pH 8.0 was used for antigen retrieval. Negative control (right) obtained from antibody was pre-absorbed by immunogen peptide.

产品名称: GRIN2 Rabbit Polyclonal Antibody
产品货号: APRab11766



Immunohistochemical analysis of paraffin-embedded Human colon cancer. Antibody was diluted at 1:100 (4°,overnight) . High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negative contrl (right) obtaned from antibody was pre-absorbed by immunogen peptide.

注意事项 (Note)

For research use only .