

产品名称: GRTP1 Rabbit Polyclonal Antibody
产品货号: APRab11805



产品概述 (Summary)

产品名称 (Production Name)	GRTP1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GRTP1
别名 (Alternative Names)	GRTP1; TBC1D6; Growth hormone-regulated TBC protein 1; TBC1 domain family member 6
基因 ID (Gene ID)	79774.0
蛋白 ID (SwissProt ID)	Q5TC63.The antiserum was produced against synthesized peptide derived from human GRTP1. AA range:71-120

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	35kDa

产品名称: GRTP1 Rabbit Polyclonal Antibody
产品货号: AP Rab11805




研究背景 (Background)

function: May act as a GTPase-activating protein for Rab family protein(s)., similarity: Contains 1 Rab-GAP TBC


domain., function: May act as a GTPase-activating protein for Rab family protein(s)., similarity: Contains 1 Rab-GAP TBC domain.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from 293 cells, using GRTP1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from Jurkat cells using GRTP1 antibody.

产品名称: GRTP1 Rabbit Polyclonal Antibody
产品货号: APRab11805



Western Blot analysis of various cells using GRTP1 Polyclonal Antibody

注意事项 (Note)

For research use only .